

Patient	RIGGS, KEVIN	Height	m	Physician	Dr. HANSEN, PETER	Received Date	29-Jun-2025
Sex	Male	Weight	kg	Hospital	Royal North Shore Hospital	Reviewed Date	30-Jun-2025
Year Of Birth (Age)	1945 (79)	BMI	----	City	Sydney		
		EOA needed to achieve an iEOA > 0.85 cm²/m²		Country	Australia		

Clinical History

Case: 31570835

MEDTRONIC ANALYSIS

ANNULUS

Diameter (mm)	22.5	x	23.6	,	23.0
	Min		Max		Mean
Perimeter (mm)	72.4		Derived Ø (mm)		23.0
Area (mm²)	415.6		Derived Ø (mm)		23.0

LVOT


Diameter (mm)	22.6	x	28.1	,	25.3
	Min		Max		Mean
Perimeter (mm)	80.6		Derived Ø (mm)		25.7
Area (mm²)	492.5		Derived Ø (mm)		25.0

VIEWS

Cusp Overlap View RAO: 8°, Caudal: 15°

3 Cusp Coplanar View


Near Cusp Overlap View



Max Ascending Aorta Diameter (mm)	31.8
Sinotubular Junction Diameter (mm)	28.1 x 28.7
	Min Max
Sinus of Valsalva Diameter (mm)	36.1 30.3 31.2
	LCC RCC NCC
Sinus of Valsalva Height (mm)	20.1 19.1 23.4
	LCC RCC NCC
Coronary Ostia Height (mm)	7.0 12.0
	Left Right

RIGHT

CIA Min Diameter (mm)	7.2	x	7.9
EIA Min Diameter (mm)	3.7	x	6.4
Femoral Min Diameter (mm)	6.5	x	6.7




LEFT

CIA Min Diameter (mm)	7.6	x	7.8
EIA Min Diameter (mm)	6.8	x	8.0
Femoral Min Diameter (mm)	5.9	x	7.2

RIGHT

Subclavian Min Diameter (mm)	----	x	----
Annular Angulation	44°		



LEFT

Subclavian Min Diameter (mm)	----	x	----
------------------------------	------	---	------

Please review images for direct aortic evaluation.

Calcium: Mild ☐ Moderate ☐ Severe ☐

VIV ADDITIONAL MEASUREMENTS

Valve to Coronary Distance (mm)	5.5	6.3
	To LCA	To RCA
Valve to STJ Distance (mm)	0.6	2.4
	LCC	RCC

Procedural Considerations

Site reports patient has a 25 mm Perimount surgical aortic valve (SAV) - Patient appears to have a 25 mm Magna Ease 3300 surgical aortic valve (SAV) - Recommend obtaining Op report for specific SAV information

Heavy MAC


Ca++ seen in LVOT under RCC, LCC and NCC

Ca++ seen in proximal RCA and LCA


Possible short segment dissection seen in REI with regions of narrowing < 5 mm

Aorta


ANNULUS




SOV DIAMETER



LVOT





STJ






Ca++ seen in LVOT under RCC, LCC and NCC

Aorta



SINUS HEIGHT




Ca++ seen in proximal LCA

Ca++ seen in proximal RCA

AORTIC ROOT




IMPLANTER'S VIEW




Cusp Overlap View






Annular Angulation



Femoral Access - Right




Femoral Access - Left





Additional Femoral Images



RAO: 0°
Caudal: 0°



Additional Images






Title


Magna Ease™ 3300 Valve¹⁰

Surgical Valve Size	19 mm	21 mm	23 mm	25 mm	27 mm	29 mm	31 mm	
Inner Diameter [†]	18	20	22	24	26	28	-	
Height	13	14	15	16	17	18	-	


Site reports patient has a 25 mm Perimount surgical aortic valve (SAV) - Patient appears to have a 25 mm Magna Ease 3300 surgical aortic valve (SAV) - Recommend obtaining Op report for specific SAV information




LCA - VTA
25 mm Magna Ease = 24 mm manufacturer's inner diameter




RCA - VTC
25 mm Magna Ease = 24 mm manufacturer's inner diameter




LCC - VTSTJ
25 mm Magna Ease = 24 mm manufacturer's inner diameter



RCC - VTSTJ
25 mm Magna Ease = 24 mm manufacturer's inner diameter



LCC/RCC - VTA at top of SAV stent posts (3 posts seen)
25 mm Magna Ease = 24 mm manufacturer's inner diameter





Patient valve selection criteria

Evolut FX bioprosthesis valve size selection

Size		23 mm	26 mm	29 mm	34 mm
Annulus diameter (A)	23.0 mm	18–20 mm	20–23 mm	23–26 mm	26–30 mm
Annulus perimeter [†]	72.4 mm	56.5–62.8 mm	62.8–72.3 mm	72.3–81.7 mm	81.7–94.2 mm
Sinus of Valsalva diameter (mean) (B)	32.5 mm	≥ 25 mm	≥ 27 mm	≥ 29 mm	≥ 31 mm
Sinus of Valsalva height (mean) (C)	20.9 mm	≥ 15 mm	≥ 15 mm	≥ 15 mm	≥ 16 mm
Oversizing Percentage		0%	13%	26%	48%

[†]Annulus perimeter = annulus diameter x π .

Selection criteria

Access consideration by MSCT

Minimum transarterial access vessel diameter

Aortic root angulation, femoral access

Aortic root angulation, left subclavian

Aortic root angulation, right subclavian

Vascular access location, direct aortic access

IFU guidance by MSCT

Evolut FX 23/26/29 mm TAVs ≥ 5.0 mm

Evolut FX 34 mm TAV ≥ 6.0 mm

Not recommended if > 70 degrees.

Not recommended if > 70 degrees.[§]

Not recommended if > 30 degrees.[§]

Ascending aorta access site ≥ 60 mm from basal plane.[¶]

[§]Patients with a patent LIMA or RIMA graft must present with access vessel diameters that are either ≥ 5.5 mm when using model D-EVOLUTFX-2329 or ≥ 6.5 mm when using model D-EVOLUTFX-34.

[¶]For direct aortic access, ensure access site and trajectory are free of patent RIMA or preexisting patent RIMA graft.

Note the position of any SVGs


(A) Annulus diameter

(B) Sinus of Valsalva diameter

(C) Sinus of Valsalva height

(D) Frame height (≈ 45 mm, not including paddles)

Illustration not to scale.



CAUTION: For distribution only in markets where CoreValve™ Evolut™ R, CoreValve™ Evolut™ PRO, Evolut™ PRO+, Evolut™ FX Systems and Evolut™ FX+ Systems are approved. See the CoreValve™ Evolut™ R, the CoreValve™ Evolut™ PRO, the Evolut™ PRO+, the Evolut™ FX and the Evolut™ FX+ device manuals for detailed information regarding the instructions for use, the implant procedure, indications, contraindications, warnings, precautions, and potential adverse events. For further information, contact your local Medtronic representative and/or consult the Medtronic website at medtronic.eu. For applicable products, consult instructions for use on manuals.medtronic.com. Manuals can be viewed using a current version of any major internet browser. For best results, use Adobe Acrobat® Reader with the browser. The commercial name of the Evolut™ R device is Medtronic CoreValve™ Evolut™ R System, the commercial name of the Evolut™ PRO device is Medtronic CoreValve™ Evolut™ PRO System, the commercial name of the Evolut™ PRO+ device is Medtronic CoreValve™ Evolut™ PRO+ System, the commercial name of the Evolut™ FX device is Medtronic Evolut™ FX System and the commercial name of the Evolut™ FX+ device is Medtronic Evolut™ FX+ System.

CAUTION: This report is provided based on information and images provided by the physician to Medtronic. This report is intended to be a resource to support physicians in their determination of proper case selection, device sizing and procedure planning, and is in no way intended to constitute medical advice or in any way replace the independent medical judgment of a trained and licensed physician with respect to any patient needs or circumstances. Physicians must conduct their own measurements and make their own medical judgments based on all of their patient's clinical and diagnostic records and images. Physician is solely responsible for all decisions and any medical judgments relating to patient diagnosis and treatment, including case selection and sizing of the device. Please see the complete Instructions of Use for all product indications, contraindications, precautions, warnings, and adverse events.

See the CoreValve™ Evolut™ R, the CoreValve™ Evolut™ PRO, the Evolut™ PRO+ device and the Evolut™ FX device manuals for detailed information regarding the instructions for use, the implant procedure, indications, contraindications, warnings, precautions, and potential adverse events. For further information, contact your local Medtronic representative and/or consult the Medtronic website at medtronic.eu.

For applicable products, consult instructions for use on manuals.medtronic.com. Manuals can be viewed using a current version of any major internet browser. For best results, use Adobe Acrobat® Reader with the browser.

The commercial name of the Evolut™ R device is Medtronic CoreValve™ Evolut™ R System, the commercial name of the Evolut™ PRO device is Medtronic CoreValve™ Evolut™ PRO System, the commercial name of the Evolut™ PRO+ device is Medtronic Evolut™ PRO+ System, and the commercial name of the Evolut™ FX device is Medtronic Evolut™ FX System.

CE
0344

Medtronic

Europe

Medtronic International Trading Sàrl.
Route du Molliau 31
Case postale
CH-1131 Tolochenaz
www.medtronic.eu

UC202202685a-evolut-fx-patient-evaluation-criteria-en-we-8664725

©2023 Medtronic. All rights reserved.

Medtronic, Medtronic logo, and Engineering the extraordinary are trademarks of Medtronic.

™*Third-party brands are trademarks of their respective owners.

All other brands are trademarks of a Medtronic company.

medtronic.eu

Notes:

Conclusion:
Reviewer Name: #42
Review Date: 30-Jun-2025